





BIOMATE & BIOMATE PLUS
DENTAL IMPLANT SYSTEM



Biomate SWISS
Implant Technology







DENTAL IMPLANT SYSTEM

Index

Product Features

About Biomate SWISS.....	1
Global Distributor.....	3
PDL [®] Surface Treatment.....	5
A Total Solution.....	7
IRB Clinical Trial	8
Biomate Implant System.....	9
Biomate Plus Implant System.....	10
Biomate Implant Design.....	11
Biomate Plus Implant Design.....	13
Biomate Dimension Table.....	15

Product Configurations

Solid Abutment / Simple Abutment

Abutment Level Impression.....	21
--------------------------------	----

Simple / Angled / Shaping / UCLA

Fixture Level Impression.....	31
-------------------------------	----

CAD CAM / Ti-base / Temporary / Premilled

Fixture Level Impression.....	37
-------------------------------	----

Multi-Unit Straight / Multi-Unit Angled

Abutment Level Impression.....	42
--------------------------------	----

Ball Abutment / Positioner

Abutment Level Impression.....	49
--------------------------------	----

Surgical Kits.....	55
--------------------	----

Surgical Instruments.....	73
---------------------------	----

Depth Marks on Biomate Instruments.....	74
---	----

Surgical Procedure.....	81
-------------------------	----

Preoperative Evaluation.....	82
------------------------------	----

Biomate Dental Implant System

The spirit of the craftsman is deeply characterized
In the soul of Biomate Swiss Implant System



Switzerland; The worldwide-recognized authority on mechanical design and craftsmanship. The craftsmanship is not only presents high precision, but also combines professional knowledge and ultimate beauty format.

The credit of Biomate SWISS implant system is because of its excellent surface treatment. This surface treatment is a collection of historical quenching, vision and view point of worldwide.



PDL® Laser Surface Treatment



Hemocompatibility



Structure of Micro-Nino Pore



Better Biocompatibility

With the primary mission of "being the best support team of dentists, and allowing patients to regain confidence and charming smiles", Biomate Swiss implant system listens closely to users' needs, provides localized and considerate services, and creates customized products and services.

We strive to become the physician's most trustworthy partner to help restore the smiles on the faces of our patients with our innovations and services.

Combining exclusive technology and professional certification, to develop PDL® (Precision Dimension Laser) core laser technology, market the Biomate SWISS implant system globally, and create a global operation model that integrates hardware and teaching services.

Global Distributor



The first medical device manufacturer in Taiwan that won a number of European innovation and invention awards.



European Exhibition of Creativity and Innovation / Special Award



Archimedes Moscow International Salon of Inventions and Innovation Technologies / Gold Medal



International Warsaw Invention Show / Silver Medal



International Salon Of Inventions and Innovation Technologies / Gold Medal



| PDL® Surface Treatment

BIOMATE Dental Implant is Designed with New Concept and New Technique, It Possesses the Best Stability to Ensure Long-Term Efficacy.

PDL® (Precision Dimension Laser) Surface Treatment applies precise parametric design. Through this high efficacy laser luminous energy, the structure can be formed on the implant surface with qualitative micro and complex 3D pores, it even carves out multiple micro-channel in linear arrangement which is suitable for the adhesion and growth of osteocytes that help increase the contact surface area of the bone and fixture, optimizing the effect of cell proliferation and osseointegration.

Better Biocompatibility

With PDL® surface treatment, a micro composite structure multiple micro-channel on the surface of the implant will guide the predecessor cells of osteoblast moving to the surface of the implant and distributing according to the track structure. The feature can increase rapid distribution and stability of the cells. When the cells are moving nearby complex 3D pores, they will randomly attach to the pores and differentiate to osteoblast which can accelerate the growth of new bone and shorten the time of osseointegration.


Completely Cleanly production process

The production interface adopts the environmental friendly manufacturing process without any chemical media. Thus, the risk of chemical material remains can be prevented.

Osseointegration


The complex micro 3D pores of BIOMATE fixture surface can effectively help to:

- Optimize the adhesion and growth of osteocytes
- Accelerate the healing of wound
- Improve osseointegration




Contact Area

PDL®(Precision Dimension Laser) Surface Treatment applies precise parametric design and through high efficacy laser luminous energy to strike qualitative micro , complex 3D texture that help increase the contact surface area of the bone and fixture, optimizes the effect of cell proliferation and osseointegration.



Cell Adhesion

The structure surface of multiple micro channel created by PDL® technique can help the adhesion of hydrophilic protein like cells, fibrin and so on. There is special metal solution and molecular arrangement in the micro pore. When osteocyte enters the pore, it can stably adhere to the structure, plus with the special pore size of Biomate which can accelerate the cells extension and differential, it can significantly improve the osseointegration (the red arrow in the picture shows the cell. SEM report shows the PDL® treated surface is highly suitable for the growth of osteocytes).




Hemocompatibility

Proven by experiment, the complex micro texture of BIOMATE-PLUS fixture surface has the best hemocompatibility, which absorbs blood rapidly to the fixture surface during implantation. This feature can effectively accelerate bone regeneration and osseointegration.



Cleanliness

Unlike SLA surface treatment that risk of sand & acidic residue after treatment, PDL® surface treatment modifies the fixture surface with laser that does not leave any toxic residue.



Surface Elements:Ti, N, O, C
Analysis of Chemical Bond:TiO₂
Evidence shows there is no residual on surface.
The surface is fully clean.

A Total Solution

Digital Surgical Guide & Customized Prosthetics Service

Biomate SWISS possesses leading CAD CAM(Computer-aided design / Computer-aided manufacture) dental facilities. This technique integrate 3D photographing and computerized digital application, helping to quickly and accurately produce different type of dental restorations.

Allow to produce digital artificial teeth product with high added-value.

Provide Customized Abutment Service

Solid Titanium
customized abutment


Zirconia plus Titanium
base customized abutment

BIOnavi. Digital Implant System

Biomate ArchFixation



IRB Clinical Trial



Description:

As shown on the above bar chart, Biomate SWISS fixture with the patent functional laser surface treatment shows best hemocompatibility.



Source of Origin:
Surface Functionalization Treatment on Dental Implant; P118-131,
2013 Clinical Oral Implantology II, Taipei Congress of Oral
Implantologists

Biomate Implant System

Suitable for:

- ✓ Surgery needs initial stability
- ✓ Immediate replacement and immediate loading
- ✓ D3~D4 bone

















Material:

- Medical grade 4 pure titanium
- In the posterior area, please use implants of over diameter 4.1mm (inclusive) size



Material: Medical grade 4 pure titanium

Unit:mm , Scale 1 : 1.5 / mm

	L8	L10	L12	L14
SD D3.3 Hex 2.0				
	1AA-001	1AA-002	1AA-003	1AA-004
SD D4.1 Hex 2.0				
	1AA-005	1AA-006	1AA-007	1AA-008
RD D4.8 Hex 2.5				
	1AA-009	1AA-010	1AA-011	1AA-012
RD D5.5 Hex 2.5				
	1AA-013	1AA-014	1AA-015	1AA-016

Biomate Plus implant System

Suitable for:

- ✓ All Bone densities(D1~D4)
- ✓ Minimally invasive surgery













Material:

- Medical grade 4 pure titanium
- In the posterior area, please use implants of over diameter 4.0mm (inclusive) size




Material: Medical grade 4 pure titanium

Unit:mm , Scale 1 : 1.5 / mm

	L8	L10	L12	L14
SD D3.5 Hex 2.0				
	1AA-017	1AA-018	1AA-019	1AA-020
SD D4.0 Hex 2.0				
	1AA-021	1AA-022	1AA-023	1AA-024
SD D4.5 Hex 2.0				
	1AA-025	1AA-026	1AA-027	1AA-028

Biomate Implant Design



Cross Section Design

- Cross section design for preventing denture rotation.




0.3mm Platform Switch Design

- The anti-bacterial and polishing design on the platform avoid the growth of dental plaque and decrease bone absorption.
- Avoid bone loss and gum recession.



Minor External Expansion Design (Ø4.1 · Ø4.8 · Ø5.5)

- Increase fixture stability in primary stage and helpful for stability in extracted socket.



Root Form Design

- Tapered body for use in anatomically constricted area.



Self-Tapping Thread Design

- Ensure better primary stability and avoid excessive force causing cortical bone absorption and promotes bone tapering during insertion.



Arced Root Design

- Arc shape at the bottom of the fixture.
- Avoid damaging vital structure like inferior alveolar nerve or maxillary sinus.



Anodized Coloring Treatment

Medical coloring treatment technique

- An oxide layer formed by anodic treatment color the abutment gold, increasing the aesthetic of the gum.
- Coloring treatment can improve the affinity of the abutment and gum, thus solidify the connection of fibro-tissue.



Various Abutment Specifications

- Various abutment configurations to suit different oral conditions.
- M1.6 abutment screw is designed to be compatible to all abutments.
- Screw Driver Hex 1.25mm is designed to match all abutment configuration (Ball abutment / Positioner abutment / Multi Unit abutment excluded).




Internal Hexagon & 10° Morse Taper

Internal Hex Design

- Avoid abutment rotation, increasing stability of the connection between the fixture and the abutment.

Abutment Screw

- The abutment screw connects the abutment firmly to the fixture; such component is not subject to load, eliminating the risk of breakage.




10° Morse Tapper Design

- Ensure firm connection between the fixture and the abutment, eliminating the possibility of unscrewing and micromovement, avoiding mechanical stimulation to surrounding tissue as well as preventing the intrusion of cell and bacteria.

Trapezoid Thread Design


- The lower part of the thread carries a larger angle of 25 degree allowing easier insertion of implant.
- The upper part carries a smaller angle of 15 degree preventing implant from dislodging.
- The space between the threads also forms an asymmetric trapezoid shape with a 0.3mm depth and 0.8mm apart. During the insertion of implant, the lower slope of this asymmetric trapezoid space (the upper part of the thread) squeezing the bone upward and compacting the bone into the space.
- The asymmetric trapezoid also benefits the laser processing on implant surface which providing good environment for optimal bone



Biomate Plus Implant Design

Single pore dimension of abutment design

- With the simplification of the restoration components for Biomate & Biomate-Plus implant system. We can complete the locked function of restoration components by only using 1.25mm Hex Driver HP/RT.



Cross Section Design

- Cross section design for preventing denture rotation.

0.3mm Platform Switch Design

- The anti-bacterial and polishing design on the platform avoid the growth of dental plaque and decrease bone absorption.
- Avoid bone loss and gum recession.

0.7mm vertical machined surface

- May adjust different insertion depth according to the requirement of implant area.

0.8mm non-continuous parallel thread

- Excellent effect for maintaining the height of bone level.

Root Form Design

- Tapered body for use in anatomically constricted area.



Self-Tapping Thread Design

- Ensure better primary stability and avoid excessive force causing cortical bone absorption and promotes bone tapering during insertion.

15°

25°

Arced Root Design

- Arc shape at the bottom of the fixture.
- Avoid damaging vital structure like inferior alveolar nerve or maxillary sinus.



0.3mm
0.7mm

Anodized Coloring Treatment

Medical coloring treatment technique

- An oxide layer formed by anodic treatment color the abutment gold, increasing the aesthetic of the gum.
- Coloring treatment can improve the affinity of the abutment and gum, thus solidify the connection of fibro-tissue.



Various Abutment Specifications

- Various abutment configurations to suit different oral conditions.
- M1.6 abutment screw is designed to be compatible to all abutments.
- Screw Driver Hex 1.25mm is designed to match all abutment configuration (Ball abutment / Positioner abutment / Multi Unit abutment excluded).




Internal Hexagon & 10° Morse Taper

Internal Hex Design

- Avoid abutment rotation, increasing stability of the connection between the fixture and the abutment.

Abutment Screw

- The abutment screw connects the abutment firmly to the fixture; such component is not subject to load, eliminating the risk of breakage.




10° Morse Tapper Design

- Ensure firm connection between the fixture and the abutment, eliminating the possibility of unscrewing and micromovement, avoiding mechanical stimulation to surrounding tissue as well as preventing the intrusion of cell and bacteria.

Trapezoid Thread Design

- The lower part of the thread carries a larger angle of 25 degree allowing easier insertion of implant.
- The upper part carries a smaller angle of 15 degree preventing implant from dislodging.
- The space between the threads also forms an asymmetric trapezoid shape with a 0.3mm depth and 0.8mm apart. During the insertion of implant, the lower slope of this asymmetric trapezoid space (the upper part of the thread) squeezing the bone upward and compacting the bone into the space.
- The asymmetric trapezoid also benefits the laser processing on implant surface which providing good environment for optimal bone








Biomate Dimension Table

Diameter 3.3 is designed for anterior area, not recommended for posterior area.

Unit : mm

Fixture						
	SD			RD		
	Diameter	3.3	4.1	4.8	5.5	
	Length	8	10	12	14	
	Platform	3.3	4.1	4.8	5.5	
	Body Diameter	2.8	2.8	3.3	4.0	
	Bevel Height	0.3	0.4	0.4	0.4	
	Final Drill	2.8 (Blue)	2.8 (Blue)	3.3 (Purple)	4.0 (Orange)	
Counter Sink	/	4.1 (Green)	4.8 (Purple)	5.5 (Orange)		
Healing Abutment						
	SD			RD		
	Diameter	4.0	4.5	5.0	5.0	6.0
	Height		2	3	5	7
Simple Abutment (Hex / Non Hex)						
	SD			RD		
	Diameter	4.0	4.5	5.0	5.0	6.0
	Height		4.0	5.5	7.0	
Gingival Height						
	1			1		
	2			2		
	3			3		
4			4			

Biomate Plus Dimension Table

Unit : mm			
Fixture	 SD	 SD	 SD
Diameter	3.5	4.0	4.5
Length	8 10 12 14	8 10 12 14	8 10 12 14
Platform	3.5 (Dark blue)	4.0 (Pink)	4.5 (Brown)
Body Diameter	2.8 (Blue)	3.3 (Purple)	4.0 (Orange)
Bevel Height	1.0	1.0	1.0
Final Drill	2.8 (Blue)	3.3 (Purple)	4.0 (Orange)
Profile Drill	3.5 (Dark blue)	4.0 (Pink)	4.5 (Brown)
Healing Abutment	 SD		
Diameter	4.0	4.5	5.0
Height	2 3 5 7		
Simple Abutment (Hex / Non Hex)	 SD		
Diameter	4.0	4.5	5.0
Height	4.0 5.5 7.0		
Gingival Height	1 2 3 4		

Product Configurations

Healing Abutment (mark)	19
Cover Screw.....	20
Membrane Screw.....	20

PROSTHETIC FLOW DIAGRAM 1

Solid Abutment / Simple Abutment Abutment Level Impression.....	21
Solid Abutment.....	22
Protect Cap.....	23
Impression Coping.....	23
Cylinder-Single/ Cylinder-Bridge.....	24
Abutment Analog.....	24
Simple Abutment.....	25
Laboratory Screw.....	27
Implant Analog.....	27
Try-in Abutment.....	28
Impression Post-Open Tray.....	29
Impression Post-Open Tray Screw.....	29
Impression Post-Close Tray.....	30
Impression Post-Close Tray Screw.....	30

PROSTHETIC FLOW DIAGRAM 2

Simple / Angled / Shaping / UCLA Fixture Level Impression.....	31
15° Angled Abutment.....	32
25° Angled Abutment.....	33
15° Angled Try-in Abutment.....	34
25° Angled Try-in Abutment.....	34
Shaping Abutment.....	35
UCLA Abutment.....	36

Product Configurations

PROSTHETIC FLOW DIAGRAM 3

CAD CAM / Ti-Base / Temporary / Premilled Fixture Level Impression.....	37
ScanBody.....	38
Ti-Base Abutment.....	38
Premilled Abutment.....	39
Temporary Abutment.....	40
Temporary Abutment (PEEK)	41


PROSTHETIC FLOW DIAGRAM 4

Multi-Unit Straight / Multi-Unit Angled Abutment Level Impression.....	42
Multi-Unit Straight Abutment.....	43
Multi-Unit 17° Angled Abutment.....	44
Multi-Unit 30° Angled Abutment.....	45
Multi-Unit Comfort Cap.....	45
Multi-Unit Titanium Cylinder.....	46
Multi-Unit Metal Cylinder.....	46
Multi-Unit Burn-out Cylinder.....	47
Multi-Unit Polishing Protector.....	47
Multi-Unit Analog.....	47
Multi-Unit Impression Coping Pick-up.....	48
Multi-Unit Impression Coping Transfer.....	48
Multi-Unit Abutment Adapter.....	48

PROSTHETIC FLOW DIAGRAM 5

Ball Abutment / Positioner Abutment Level Impression.....	49
Ball Abutment.....	50
Housing Retainer with O-ring / Housing with O-ring / Ball Abutment Analog / Ball Abutment Driver	51
Positioner Abutment / Positioner Torque Driver / Positioner Core Tool	52
Positioner Male Processing Kit.....	53
Positioner Replacement Male.....	53
Positioner Extended Replacement Male.....	53
Positioner Black Processing Male.....	54
Positioner Block Out Spacers.....	54
Positioner Impression Coping.....	54
Positioner Lab Analog	54


Healing Abutment (mark)



- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:10Ncm

Material:

- Medical Grade 4 Pure Titanium



D \ H	2.0	3.0	5.0	7.0
4.0	1AA-316	1AA-313	1AA-314	1AA-315
4.5	1AA-317	1AA-301	1AA-302	1AA-303
5.0	1AA-318	1AA-307	1AA-308	1AA-309
6.0	1AA-324	1AA-321	1AA-322	1AA-323



D \ H	2.0	3.0	5.0	7.0
5.0	1AA-319	1AA-304	1AA-305	1AA-306
6.0	1AA-320	1AA-310	1AA-311	1AA-312

Cover Screw

✓ Hand tightened with 1.25 hex driver

- Material:
- Medical Grade 4 Pure Titanium



Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5

SD 1AA-101



Implant Size
Ø4.8 Ø5.5

RD 1AA-102

Membrane Screw

- ✓ Used for securing membranes to implant.
✓ Threaded into the cover screw inner thread.
✓ Hand tightened with 1.25 hex driver

- Material:
- Medical Grade 4 Pure Titanium



Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5

D	G/H	0	0.5	1	2
	H	5.5	6.0	6.5	7.5
3.3		1AA-111	1AA-112	1AA-113	1AA-114
5.0		1AA-115	1AA-116	1AA-117	1AA-118




Implant Size
Ø4.8 Ø5.5


D	G/H	0	0.5	1	2
	H	6.2	6.7	7.2	8.2
4.8		1AA-103	1AA-104	1AA-105	1AA-106
6.0		1AA-107	1AA-108	1AA-109	1AA-110

Solid Abutment / Simple Abutment

Abutment Level Impression



Solid Abutment







- ✓ Abutment for producing cement-retained/combination prosthesis
- ✓ Coating:SD/Semi-Golden; RD/Golden
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:25Ncm

Material:

- Medical Grade 4 Pure Titanium





S

Implant Size
 Ø3.3 Ø4.1
 Ø3.5 Ø4.0 Ø4.5

D \ G/H		1	2	3	4
H					
Non Hex 4.0	4.0	4AA-A41	4AA-A44	4AA-A47	4AA-A50
Non Hex 4.0	5.5	4AA-A42	4AA-A45	4AA-A48	4AA-A51
Non Hex 4.0	7.0	4AA-A43	4AA-A46	4AA-A49	4AA-A52
Non Hex 4.5	4.0	4AA-A53	4AA-A55	4AA-A57	4AA-A59
Non Hex 4.5	5.5	4AA-A01	4AA-A02	4AA-A03	4AA-A04
Non Hex 4.5	7.0	4AA-A54	4AA-A56	4AA-A58	4AA-A60
Non Hex 5.0	4.0	4AA-A05	4AA-A08	4AA-A11	4AA-A14
Non Hex 5.0	5.5	4AA-A06	4AA-A09	4AA-A12	4AA-A15
Non Hex 5.0	7.0	4AA-A07	4AA-A10	4AA-A13	4AA-A16

R

Implant Size
 Ø4.8 Ø5.5

D \ G/H		1	2	3	4
H					
Non Hex 5.0	4.0	4AA-A17	4AA-A20	4AA-A23	4AA-A26
Non Hex 5.0	5.5	4AA-A18	4AA-A21	4AA-A24	4AA-A27
Non Hex 5.0	7.0	4AA-A19	4AA-A22	4AA-A25	4AA-A28
Non Hex 6.0	4.0	4AA-A29	4AA-A32	4AA-A35	4AA-A38
Non Hex 6.0	5.5	4AA-A30	4AA-A33	4AA-A36	4AA-A39
Non Hex 6.0	7.0	4AA-A31	4AA-A34	4AA-A37	4AA-A40

Solid Abutment Components

Protect Cap

- ✓ Used for Solid/Simple abutment protection and reducing patient discomfort
- ✓ Used as a temporary crown base



Implant Size

Ø3.3	Ø4.1	
Ø3.5	Ø4.0	Ø4.5
Ø4.8	Ø5.5	

D \ H	4.0	5.5	7.0
4.0	6AA-087	6AA-051	6AA-049
4.5	6AA-052	6AA-019	6AA-050
5.0	6AA-053	6AA-020	6AA-021
6.0	6AA-088	6AA-022	6AA-023

Impression Coping

- ✓ Components for Solid/Simple abutment impression
- ✓ The top lug is designed to align with the cross section of the body of abutment for accurate positioning



Implant Size




Ø3.3	Ø4.1	
Ø3.5	Ø4.0	Ø4.5
Ø4.8	Ø5.5	

D \ H	10.0
4.0	6AA-048
4.5	6AA-016
5.0	6AA-017
6.0	6AA-018

Solid Abutment Components

Cylinder-Single / Cylinder-Bridge

- ✓ Enabling the production of coping with abutment analog
- ✓ Used after casting, after cleaning the margin for proper fitting

		Cylinder-Single		Cylinder-Bridge
	D \ H	10.0	D \ H	10.0
 Implant Size	4.5	6AA-010	4.5	6AA-013
	5.0	6AA-011	5.0	6AA-014
	6.0	6AA-012	6.0	6AA-015

Ø3.3

Ø3.5

Ø4.8

Ø4.1



Ø4.0

Ø5.5

Ø4.5

Abutment Analog

- ✓ Solid/Simple abutment reproduction on model after impression

		4.0		5.5		7.0	
D \ H		4.0		5.5		7.0	
 Implant Size	4.0	6AA-054		6AA-055		6AA-056	
	4.5	6AA-057		6AA-003		6AA-058	
	5.0	6AA-004		6AA-005		6AA-006	
	6.0	6AA-007		6AA-008		6AA-009	

Ø3.3

Ø3.5

Ø4.8


Ø4.1

Ø4.0

Ø5.5

Ø4.5

Simple Abutment



Simple abutment is the same as solid abutment above margin, so all components of solid abutment can be shared with simple abutment

- ✓ Abutment for producing cement-retained/combination prosthesis
- ✓ Coating:SD/Semi-Golden; RD/Golden
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:30Ncm

Material:
• Medical Grade 4 Pure Titanium

S

Implant Size





Ø3.3





Ø3.5

Ø4.1

Ø4.0

Ø4.5

D		G/H	1	2	3	4
		H				
Hex	4.0	4.0	4AA-B41	4AA-B44	4AA-B47	4AA-B50
Hex	4.0	5.5	4AA-B42	4AA-B45	4AA-B48	4AA-B51
Hex	4.0	7.0	4AA-B43	4AA-B46	4AA-B49	4AA-B52
Hex	4.5	4.0	4AA-B53	4AA-B55	4AA-B57	4AA-B59
Hex	4.5	5.5	4AA-B01	4AA-B02	4AA-B03	4AA-B04
Hex	4.5	7.0	4AA-B54	4AA-B56	4AA-B58	4AA-B60
Hex	5.0	4.0	4AA-B05	4AA-B08	4AA-B11	4AA-B14
Hex	5.0	5.5	4AA-B06	4AA-B09	4AA-B12	4AA-B15
Hex	5.0	7.0	4AA-B07	4AA-B10	4AA-B13	4AA-B16

D		G/H	1	2	3	4
		H				
Non Hex	4.0	4.0	4AA-C41	4AA-C44	4AA-C47	4AA-C50
Non Hex	4.0	5.5	4AA-C42	4AA-C45	4AA-C48	4AA-C51
Non Hex	4.0	7.0	4AA-C43	4AA-C46	4AA-C49	4AA-C52
Non Hex	4.5	4.0	4AA-C53	4AA-C55	4AA-C57	4AA-C59
Non Hex	4.5	5.5	4AA-C01	4AA-C02	4AA-C03	4AA-C04
Non Hex	4.5	7.0	4AA-C54	4AA-C56	4AA-C58	4AA-C60
Non Hex	5.0	4.0	4AA-C05	4AA-C08	4AA-C11	4AA-C14
Non Hex	5.0	5.5	4AA-C06	4AA-C09	4AA-C12	4AA-C15
Non Hex	5.0	7.0	4AA-C07	4AA-C10	4AA-C13	4AA-C16

S

Implant Size

Ø3.3


Ø3.5

Ø4.1

Ø4.0

Ø4.5

Simple Abutment







Simple abutment is the same as solid abutment above margin, so all components of solid abutment can be shared with simple abutment

- ✓ Abutment for producing cement-retained/combination prosthesis
- ✓ Coating:SD/Semi-Golden; RD/Golden
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:30Ncm





Material:

- Medical Grade 4 Pure Titanium



D		G/H	1	2	3	4
		H				
Hex	5.0	4.0	4AA-B17	4AA-B20	4AA-B23	4AA-B26
Hex	5.0	5.5	4AA-B18	4AA-B21	4AA-B24	4AA-B27
Hex	5.0	7.0	4AA-B19	4AA-B22	4AA-B25	4AA-B28
Hex	6.0	4.0	4AA-B29	4AA-B32	4AA-B35	4AA-B38
Hex	6.0	5.5	4AA-B30	4AA-B33	4AA-B36	4AA-B39
Hex	6.0	7.0	4AA-B31	4AA-B34	4AA-B37	4AA-B40



		D	G/H	1	2	3	4
		H					
							
Non Hex	5.0	4.0	4AA-C17	4AA-C20	4AA-C23	4AA-C26	
Non Hex	5.0	5.5	4AA-C18	4AA-C21	4AA-C24	4AA-C27	
Non Hex	5.0	7.0	4AA-C19	4AA-C22	4AA-C25	4AA-C28	
Non Hex	6.0	4.0	4AA-C29	4AA-C32	4AA-C35	4AA-C38	
Non Hex	6.0	5.5	4AA-C30	4AA-C33	4AA-C36	4AA-C39	
Non Hex	6.0	7.0	4AA-C31	4AA-C34	4AA-C37	4AA-C40	

Simple Abutment Components

Abutment Screw

✓ Screw used to assemble implant and abutment

- Material:
- Medical Grade 5 Titanium Alloy



Implant Size

Ø3.3	Ø4.1
Ø3.5	Ø4.0
Ø4.8	Ø5.5



Laboratory Screw

✓ Abutment screw for lab work


✓ Coating:Golden

- Material:
- Medical Grade 5 Titanium Alloy



Implant Size

Ø3.3	Ø4.1
Ø3.5	Ø4.0
Ø4.8	Ø5.5




Implant Analog

✓ Lab analog for fixture level impression



Implant Size

Ø3.3	Ø4.1
Ø3.5	Ø4.0
	Ø4.5




Implant Size

Ø4.8	Ø5.5
------	------



Simple Abutment Components

Try-in Abutment



- ✓ Used in selecting specifications of Solid/Simple abutment
- ✓ Coating:SD/Semi-Golden; RD/Golden

Material:

- Medical Grade 4 Pure Titanium



Implant

Size

Ø3.3 Ø4.1

Ø3.5 Ø4.0 Ø4.5

D

4.0

4.5

5.0

3AA-054

3AA-046

3AA-047



Implant

Size

Ø4.8 Ø5.5

D

5.0


6.0

3AA-048

3AA-049


Impression Post-Open Tray

- ✓ Components for fixture level impression taking with open tray
- ✓ Multi-undercutting design that is stably fixed within the impression body
- ✓ Hand tightened with 1.25 hex driver




Implant Size

Ø3.3	Ø4.1
Ø3.5	Ø4.0
	Ø4.5




D \ H	10.8
Hex 4.0	6AA-024
Hex 5.0	6AA-025
Non Hex 4.0	6AA-035
Non Hex 5.0	6AA-036



Implant Size


Ø4.8	Ø5.5
------	------



D \ H	10.8
Hex 5.0	6AA-059
Hex 6.0	6AA-060
Non Hex 5.0	6AA-063
Non Hex 6.0	6AA-064


Impression Post-Open Tray Screw

- ✓ Screw used to assemble implant and impression post-open tray



Implant Size

Ø3.3	Ø4.1
Ø3.5	Ø4.0
	Ø4.5
Ø4.8	Ø5.5



L	22	26
	6AA-028	6AA-029

Impression Post-Close Tray

- ✓ Components for fixture level impression taking with closed tray
- ✓ Undercutting design for stable fastening and accurate repositioning
- ✓ Hand tightened with 1.25 hex driver



Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5

D \ H		10.5
Hex	4.0	6AA-030
Hex	5.0	6AA-031
Non Hex	4.0	6AA-039
Non Hex	5.0	6AA-040



Implant Size
Ø4.8 Ø5.5

D \ H		10.5
Hex	5.0	6AA-061
Hex	6.0	6AA-062
Non Hex	5.0	6AA-065
Non Hex	6.0	6AA-066

Impression Post-Close Tray Screw

- ✓ Screw used to assemble implant and impression post-close tray




Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5
Ø4.8 Ø5.5


L		16.5
		6AA-034

Simple / Angled / Shaping / UCLA

Fixture Level Impression



15° Angled Abutment



- ✓ Used when a prosthesis's path adjustment is necessary at 15°
- ✓ Coating:SD/Semi-Golden; RD/Golden
- ✓ The angled direction is pointed to hexagonal edge
- ✓ Accurate specification selected by angled Try-in Abutment
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:30Ncm

Material:

- Medical Grade 4 Pure Titanium



Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5


D		G/H	1	2	3	4
	H					
Hex	4.0	8.5	4AA-D17	4AA-D01	4AA-D18	4AA-D02
Hex	5.0	8.5	4AA-D19	4AA-D03	4AA-D20	4AA-D04
Non Hex	4.0	8.5	4AA-D25	4AA-D09	4AA-D26	4AA-D10
Non Hex	5.0	8.5	4AA-D27	4AA-D11	4AA-D28	4AA-D12



Implant Size
Ø4.8 Ø5.5

D		G/H	1	2	3	4
	H					
Hex	5.0	8.5	4AA-D21	4AA-D05	4AA-D22	4AA-D06
Hex	6.0	8.5	4AA-D23	4AA-D07	4AA-D24	4AA-D08
Non Hex	5.0	8.5	4AA-D29	4AA-D13	4AA-D30	4AA-D14
Non Hex	6.0	8.5	4AA-D31	4AA-D15	4AA-D32	4AA-D16

25°Angled Abutment



- ✓ Used when a prosthesis's path adjustment is necessary at 25°
- ✓ Coating:SD/Semi-Golden; RD/Golden
- ✓ The angled direction is pointed to hexagonal edge
- ✓ Accurate specification selected by angled Try-in Abutment
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:30Ncm

Material:

- Medical Grade 4 Pure Titanium



Implant
Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5

D \ G/H		1	2	3	4
H					
Hex 4.0	8.5	4AA-E17	4AA-E01	4AA-E18	4AA-E02
Hex 5.0	8.5	4AA-E19	4AA-E03	4AA-E20	4AA-E04
Non Hex 4.0	8.5	4AA-E25	4AA-E09	4AA-E26	4AA-E10
Non Hex 5.0	8.5	4AA-E27	4AA-E11	4AA-E28	4AA-E12




Implant
Size
Ø4.8 Ø5.5

D \ G/H		1	2	3	4
H					
Hex 5.0	8.5	4AA-E21	4AA-E05	4AA-E22	4AA-E06
Hex 6.0	8.5	4AA-E23	4AA-E07	4AA-E24	4AA-E08
Non Hex 5.0	8.5	4AA-E29	4AA-E13	4AA-E30	4AA-E14
Non Hex 6.0	8.5	4AA-E31	4AA-E15	4AA-E32	4AA-E16

Angled Abutment Components

15° Angled Try-in Abutment



✓ Used in selecting diameter, and G/H of angled abutment when in oral cavity or in model

Material:

- Medical Grade 4 Pure Titanium



S

Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5

D \ G/H	2	4
4.0	3AA-068	3AA-072
5.0	3AA-070	3AA-071

R

Implant Size
Ø4.8 Ø5.5

D \ G/H	2	4
5.0	3AA-069	3AA-073
6.0	3AA-074	3AA-075

25° Angled Try-in Abutment



S

Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5


D \ G/H	2	4
4.0	3AA-076	3AA-077
5.0	3AA-078	3AA-079

R

Implant Size
Ø4.8 Ø5.5

D \ G/H	2	4
5.0	3AA-080	3AA-081
6.0	3AA-082	3AA-083




Shaping Abutment




- ✓ Used when an abutment's path must be altered or a prosthesis's margin area must be customized
- ✓ Coating:SD/Semi-Golden; RD/Golden
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:30Ncm

Material:




- Medical Grade 4 Pure Titanium


D \ G/H		0	1.5	3
H				
Hex 4.0	11.0	4AA-H41	4AA-H42	4AA-H43
Hex 5.0	11.0	4AA-H44	4AA-H45	4AA-H46
Non Hex 4.0	11.0	4AA-H53	4AA-H54	4AA-H55
Non Hex 5.0	11.0	4AA-H56	4AA-H57	4AA-H58



Implant Size

Ø3.3	Ø4.1
Ø3.5	Ø4.0
	Ø4.5


D \ G/H		0	1.5	3
H				
Hex 5.0	11.0	4AA-H47	4AA-H48	4AA-H49
Hex 6.0	11.0	4AA-H50	4AA-H51	4AA-H52
Non Hex 5.0	11.0	4AA-H59	4AA-H60	4AA-H61
Non Hex 6.0	11.0	4AA-H62	4AA-H63	4AA-H64



Implant Size

Ø4.8	Ø5.5
------	------

UCLA Abutment



- ✓ Used in producing cement-retained/combination/screw-retained prosthesis
- ✓ Used when path, aesthetics, or space have limitations
- ✓ After customization, prosthesis must be produced by casting using dental-quality CCM
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:30Ncm

Material:

- Upper part/POM; Lower part/ CCM(Co-Cr-Mo)



Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5



Hex	4.5	10.0	4AA-K01
Non Hex	4.5	10.0	4AA-K03




Implant Size
Ø4.8 Ø5.5



Hex	5.0	10.0	4AA-K02
Non Hex	5.0	10.0	4AA-K04


CAD CAM / Ti-Base / Temporary / Premilled

Fixture Level Impression




ScanBody

- ✓ Scan body for intra oral scan
 - ✓ Hand tightened with 1.25 hex driver
- Material:
- Medical Grade 4 Pure Titanium



Ti-Base Abutment


- ✓ Abutment for producing cement-retained/combination/ screw-retained prosthesis
 - ✓ Used for producing Ti+Zr custom abutment with CAD/CAM equipment
 - ✓ Biomet's official implant library provided
 - ✓ Use fixture level or intra oral scan body impression
 - ✓ Tightened with 1.25 hex driver
 - ✓ Recommended tightening torque:30Ncm
- Material:
- Medical Grade 4 Pure Titanium



D	GH	0.5			1.5			2.5		
		H			H			H		
4.0	4.0	4AA-L01			4AA-L03			4AA-L05		
6.0	4.0	4AA-L07			4AA-L09			4AA-L11		

D	GH	0.5			1.5			2.5		
		H			H			H		
4.0	4.0	4AA-L02			4AA-L04			4AA-L06		
6.0	4.0	4AA-L08			4AA-L10			4AA-L12		

Premilled Abutment



- ✓ Milling equipment for dental work to product custom abutment
- ✓ Biomate's offical implant library provided
- ✓ Use fixture level or intra oral scan body impression
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:30Ncm

Material:

- Medical Grade 4 Pure Titanium



Implant Size
 Ø3.3 Ø4.1
 Ø3.5 Ø4.0 Ø4.5

D \ H 20.0



10.0 4AA-P01

D \ H 20.0




10.0 4AA-P02



Implant Size
 Ø4.8 Ø5.5

Temporary Abutment





- ✓ Used in producing temporary prosthesis
- ✓ Structure enabling easy customization
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:30Ncm

Material:

- Medical Grade 4 Pure Titanium





Implant
Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5


D \ G/H		0.5	2
	H		
Hex 4.0	10.0	4AA-G05	4AA-G07
Hex 4.5	10.0	4AA-G01	4AA-G08
Non Hex 4.0	10.0	4AA-G06	4AA-G10
Non Hex 4.5	10.0	4AA-G03	4AA-G11



Implant
Size
Ø4.8 Ø5.5

D \ G/H		0.5	2
	H		
Hex 5.0	10.0	4AA-G02	4AA-G09
Non Hex 5.0	10.0	4AA-G04	4AA-G12



Temporary Abutment (PEEK)





- ✓ Used in producing temporary prosthesis
- ✓ Structure enabling easy customization
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:30Ncm

Material:


- Medical PEEK

		D		G/H	0.5	2
				H		
<div>S</div> <div>Implant Size</div> <div> <div> <div>Ø3.3</div> <div>Ø3.5</div> </div> <div> <div>Ø4.1</div> <div>Ø4.0</div> <div>Ø4.5</div> </div> </div>	Hex	4.0	10.0		4AA-G13	4AA-G19
	Hex	4.5	10.0		4AA-G14	4AA-G20
	Non Hex	4.0	10.0		4AA-G16	4AA-G22
	Non Hex	4.5	10.0		4AA-G17	4AA-G23


		D		G/H	0.5	2
				H		
<div>R</div> <div>Implant Size</div> <div> <div>Ø4.8</div> <div>Ø5.5</div> </div>	Hex	5.0	10.0		4AA-G15	4AA-G21
	Non Hex	5.0	10.0		4AA-G18	4AA-G24

Multi-Unit Straight / Multi-Unit Angled


Abutment Level Impression








Multi-Unit Straight Abutment








- ✓ Used for producing screw-retained prosthesis in multiple case
 - ✓ Tightened with multi-unit adapter
 - ✓ Coating:SD/Semi-Golden; RD/Golden
 - ✓ Recommended tightening torque:Single/30Ncm; Biomate Archfixation/35Ncm
- Material:
- Medical Grade 5 Titanium Alloy




G/H	1	2	3	4	5
Non Hex					
	4AA-M01	4AA-M02	4AA-M03	4AA-M04	4AA-M05



G/H	1	2	3	4	5
Non Hex					
	4AA-M06	4AA-M07	4AA-M08	4AA-M09	4AA-M10

Multi-Unit 17° Angled Abutment



- ✓ Used for producing screw-retained prosthesis in multiple case
- ✓ Up to 60°path compensation(two implant standard)
- ✓ Coating:SD/Semi-Golden; RD/Golden
- ✓ Recommended tightening torque:
Single/30Ncm; Biomate Archfixation/15Ncm

Material:

- Medical Grade 5 Titanium Alloy

S

Implant Size
 Ø3.3 Ø4.1
 Ø3.5 Ø4.0 Ø4.5

Hex
Non Hex

G/H	2	3	4
	4AA-T02	4AA-T03	4AA-T04
	4AA-T06	4AA-T07	4AA-T08


R

Implant Size
 Ø4.8 Ø5.5

Hex
Non Hex

G/H	2	3	4
	4AA-T10	4AA-T11	4AA-T12
	4AA-T14	4AA-T15	4AA-T16

Multi-Unit 30° Angled Abutment



- ✓ Used for producing screw-retained prosthesis in multiple case
- ✓ Up to 60°path compensation(two implant standard)
- ✓ Coating:SD/Semi-Golden; RD/Golden
- ✓ Recommended tightening torque:
Single/30Ncm; Biomate Archfixation/15Ncm

Material:

- Medical Grade 5 Titanium Alloy



Implant

Size

Ø3.3

Ø3.5

Ø4.1

Ø4.0

Ø4.5

Hex
Non Hex

G/H

3

4

5



4AA-U02

4AA-U06



4AA-U03

4AA-U07



4AA-U04

4AA-U08



Implant

Size

Ø4.8

Ø5.5

Hex
Non Hex

G/H

3

4

5



4AA-U10

4AA-U14



4AA-U11

4AA-U15



4AA-U12

4AA-U16

Multi-Unit Straight / Multi-Unit Angled Components

Multi-Unit Comfort Cap

- ✓ Used for protecting multi-unit abutment in the oral cavity
- ✓ Hand tightened with 1.25 hex driver



Implant

Size

Ø3.3

Ø3.5

Ø4.1

Ø4.0

Ø4.5

Ø4.8

Ø5.5



BSMUCC480600A

Multi-Unit Straight / Multi-Unit Angled Components

Multi-Unit Titanium Cylinder / Multi-Unit Titanium Cylinder (S)

- ✓ Used for producing temporary prosthesis in multi-unit abutment
- ✓ (S) specification suitable for overdenture with thinner diameter
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:10-15Ncm

Material:

- Medical Grade 4 Pure Titanium

S

R

Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5
Ø4.8 Ø5.5

Hex

BSMUTC48000HA

BSMUTC4800SHA

Non Hex

BSMUTC48000NA

BSMUTC4800SNA

Multi-Unit Metal Cylinder

- ✓ Used for producing screw-retained prosthesis in multi-unit abutment
- ✓ Used to produce customized prosthesis by casting with CCM
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:10-15Ncm

Material:

- Upper part/POM; Lower part/ CCM(Co-Cr-Mo)

S

R

Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5
Ø4.8 Ø5.5

Hex

BSMUMC48000HA

Non Hex

BSMUMC48000NA


Multi-Unit Straight / Multi-Unit Angled Components

Multi-Unit Burn-out Cylinder


- ✓ Used for producing screw-retained prosthesis in multi-unit abutment
- ✓ Used to produce customized prosthesis by casting with nonprecious metal alloy
- ✓ Tightened with 1.25 hex driver
- ✓ Recommended tightening torque:10-15Ncm

Material:

- POM



Hex
Non-Hex



BSMUBC48000HA
BSMUBC48000NA

Multi-Unit Polishing Protector

- ✓ Protecting the joint in the polishing procedure after producing a prosthesis using multi-unit metal/burn-out cylinder
- ✓ Hand tightened with 1.25 hex driver

Material:

- Medical Grade 4 Pure Titanium



BSMUPP480000A

Multi-Unit Analog

- ✓ Lab analog for multi-unit abutment
- ✓ Hand tightened with 1.25 hex driver




BSMUAL480000A

Multi-Unit Straight / Multi-Unit Angled Components



Multi-Unit Impression Coping Pick-up

- ✓ Components for multi-unit abutment impression taking with open tray
- ✓ Multi-undercutting design that is stably fixed within the impression body
- ✓ Hand tightened with 1.25 hex driver





Multi-Unit Impression Coping Transfer

- ✓ Components for multi-unit abutment impression taking with close tray
- ✓ Undercutting design for stable fastening and accurate repositioning
- ✓ Hand tightened with 1.25 hex driver




Multi-Unit Abutment Adapter

- ✓ Dedicated torque driver for multi-unit abutment




Ball Abutment / Positioner

Abutment Level Impression



Ball Abutment



- ✓ Abutment for overdenture with o-ring attachment
- ✓ Insertion angle compensated up to 20°
- ✓ Coating:SD/Silver; RD/Golden
- ✓ Tightened with ball abutment driver
- ✓ Recommended tightening torque:25Ncm

Material:

- Medical Grade 4 Pure Titanium



Implant

Size

Ø3.3

Ø4.1

Ø3.5

Ø4.0

Ø4.5

D	G/H	1	2	3	4	5	6
H							
Non Hex	2.25 3.3	4AA-J01	4AA-J02	4AA-J03	4AA-J04	4AA-J05	4AA-J06



Implant

Size

Ø4.8

Ø5.5

D	G/H	1	2	3	4	5	6
H							
Non Hex	2.25 3.3	4AA-J13	4AA-J14	4AA-J15	4AA-J16	4AA-J17	4AA-J18



Ball Abutment Components

Housing Retainer with O-ring ✓ Used when vertical dimension is shorter than the housing cap

S

R

Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5
Ø4.8 Ø5.5



Housing Retainer with O-ring O-ring


6AA-044 6AA-082 / 6AA-045

Housing with O-ring ✓ O-ring attachment for ball abutment
✓ O-ring replaced in metal housing

S

R

Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5
Ø4.8 Ø5.5



Housing with O-ring

6AA-067

Ball Abutment Analog ✓ Lab analog for ball abutment

S

R

Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5
Ø4.8 Ø5.5



D 2.25




6AA-043

Ball Abutment Driver ✓ Dedicated driver for ball abutment

S


R

Implant Size
Ø3.3 Ø4.1
Ø3.5 Ø4.0 Ø4.5
Ø4.8 Ø5.5



3AA-050 3AA-051 3AA-053


Positioner Abutment



- ✓ Achieves low vertical dimension, stability and various attachments with retention
- ✓ Possible path compensation up to 40°(two implant standard)
- ✓ Tightened with positioner torque driver
- ✓ Recommended tightening torque:30Ncm

Material:

- Medical Grade 5 Titanium Alloy



Positioner Abutment Components

Positioner torque driver

- ✓ Dedicated driver for positioner abutment



3AA-085

3AA-086

Positioner Core Tool

- ✓ Used in attaching and changing replacement males




3AA-087

Positioner Abutment Components


Positioner Male Processing Kit

- ✓ Component
 - Block out spacer/ denture cap connected black processing male
 - Replacement male blue/pink/clear
- ✓ Used by selecting the male with the adequate retention force for each case
- ✓ Positioner core tool for replacing the male




Positioner Replacement Male

- ✓ Retention: Approximately Blue:1.5 lbs/Pink:3 lbs/Clear:5 lbs
- ✓ Placement angled compensation up to 20°(two implant standard)
- ✓ Packing unit: single color replacement male 4ea



Positioner Extended Replacement Male


- ✓ Retention: Approximately Red:1 lbs/Orange:2 lbs/Green:4 lbs
- ✓ Placement angled compensation up to 20°~40°(two implant standard)
- ✓ Packing unit: single color replacement male 4ea



Positioner Abutment Components


Positioner Black Processing Male

- ✓ Male used in prosthesis fabrication process
- ✓ Packing unit:4ea



Positioner Block Out Spacers

- ✓ Used for sealing of the space between the abutment and the denture cap when attaching the the overdenture and denture cap in the oral cavity
- ✓ Packing unit:20ea




Positioner Impression Coping

- ✓ Pick-up impression coping for positioner abutment with close tray
- ✓ Packing unit:4ea




Positioner Lab Analog

- ✓ Lab analog for positioner abutment
- ✓ Packing unit:4ea



*Perfect craftsmanship
Safe and reliable*





Surgical Kits


BIOMATE PLUS SC KIT.....	57
BIOMATE MASTER KIT.....	59
BIOMATE FULL KIT.....	61
BIONavi MASTER KIT.....	63
BIOMATE SL KIT.....	65
BIOMATE PROSTHETIC KIT	66
Positioner / Multi-Unit KIT.....	67
Bone Expander KIT.....	68
Sinus Crestal Approach KIT	69


Surgical Kits

BIOMATE PLUS SC KIT

Color management:

- ◆ BIOMATE/ Blue(ø3.3) · Green(ø4.1) · Purple(ø4.8) · Orange(ø5.5)
- ◆ BIOMATE PLUS/ Dark Blue(ø3.5) · Pink(ø4.0) · Brown(ø4.5)





Implant Driver HP



Implant Driver RT



Screw Driver HP
Screw Driver RT

Description		Catalog No.	
1	BIOMATE PLUS SC KIT	3AA-137	Full Instruments
2	Lance Drill D1.8	3AA-038	1EA
3	Initial Drill D2.0-S	3AA-006	1EA
	Initial Drill D2.0-L	3AA-007	1EA
4	Round Bur D2.3	3AA-001	1EA
5	Paralleling Pin D2.0- D2.8	3AA-052	3EA
6	Lindemann Drill (Pilot Drill) D2.0xL 14	3AA-037	1EA
7	Drill Extender	3AA-035	1EA
8	Final Drill D2.8-S for ϕ 3.3/4.1 Implant	3AA-008	1EA
	Final Drill D2.8-L for ϕ 3.3/4.1 Implant	3AA-009	1EA
	Final Drill D3.3-S for ϕ 4.8 Implant	3AA-010	1EA
	Final Drill D3.3-L for ϕ 4.8 Implant	3AA-011	1EA
	Final Drill D4.0-S for ϕ 5.5 Implant	3AA-012	1EA
	Final Drill D4.0-L for ϕ 5.5 Implant	3AA-013	1EA
	Counter Sink D4.1	3AA-014	1EA
9	Counter Sink D4.8	3AA-015	1EA
	Counter Sink D5.5	3AA-016	1EA
10	Profile Drill D3.5	3AA-065	1EA
	Profile Drill D4.0	3AA-066	1EA
	Profile Drill D4.5	3AA-067	1EA
11	Taps D3.3 for ϕ 3.3/4.1 Implant	3AA-017	1EA
	Taps D4.0 for ϕ 4.0/4.8 Implant	3AA-018	1EA
	Taps D4.7 for ϕ 4.5/5.5 Implant	3AA-019	1EA
12	Implant Driver Hex2.0 HP-L	3AA-056	1EA
	Implant Driver Hex2.0 RT-S	3AA-030	1EA
	Implant Driver Hex2.0 RT-L	3AA-039	1EA
	Implant Driver Hex2.5 HP-L	3AA-057	1EA
	Implant Driver Hex2.5 RT-S	3AA-032	1EA
	Implant Driver Hex2.5 RT-L	3AA-040	1EA
13	Screw Driver Hex1.25-HP-L	3AA-041	1EA
	Screw Driver Hex1.25-RT-S	3AA-042	1EA
	Screw Driver Hex1.25-RT-L	3AA-043	1EA
14	Stopper L6 for Drill D2.0/D2.8	3AA-020	1EA
	Stopper L8 for Drill D2.0/D2.8	3AA-021	1EA
	Stopper L10 for Drill D2.0/D2.8	3AA-022	1EA
	Stopper L12 for Drill D2.0/D2.8	3AA-023	1EA
	Stopper L14 for Drill D2.0/D2.8	3AA-024	1EA
	Stopper L8 for Drill D3.3/D4.0	3AA-025	1EA
	Stopper L10 for Drill D3.3/D4.0	3AA-026	1EA
	Stopper L12 for Drill D3.3/D4.0	3AA-027	1EA
15	Stopper L14 for Drill D3.3/D4.0	3AA-028	1EA
	Torque Ratchet 10-40Ncm	3AA-034	1EA


BIOMATE MASTER KIT

Color management:

- ◆ BIOMATE/ Blue(ø3.3) · Green(ø4.1) · Purple(ø4.8) · Orange(ø5.5)
- ◆ BIOMATE PLUS/ Dark Blue(ø3.5) · Pink(ø4.0) · Brown(ø4.5)



Round Bur Lance Drill Initial Drill Final Drill



Lindemann Drill
(Pilot Drill)


Drill Extender

Paralleling Pin

Counter Sink D4.1
Profile Drill D3.5

Counter Sink D4.8
Profile Drill D4.0

Counter Sink D5.5
Profile Drill D4.5



Implant Driver HP



Implant Driver RT



Screw Driver HP
Screw Driver RT



Taps D3.3 for Ø3.3/4.1 Implant
Taps D4.0 for Ø4.0/4.8 Implant
Taps D4.7 for Ø4.5/5.5 Implant








Torque Ratchet

	Description	Catalog No.	
1	BIOMATE MASTER KIT	3AA-154	Full Instruments
2	Lance Drill D1.8	3AA-038	1EA
	Initial Drill D2.0-8mm	3AA-184	1EA
	Initial Drill D2.0-10mm	3AA-185	1EA
3	Initial Drill D2.0-12mm	3AA-186	1EA
	Initial Drill D2.0-14mm	3AA-187	1EA
4	Round Bur D2.3	3AA-001	1EA
5	Paralleling Pin D2.0-D2.8	3AA-052	2EA
6	Lindemann Drill (Pilot Drill) D2.0xL14	3AA-037	1EA
7	Drill Extender	3AA-035	1EA
	Final Drill D2.8-8mm for Ø3.3/4.1 Implant	3AA-188	1EA
	Final Drill D2.8-10mm for Ø3.3/4.1 Implant	3AA-189	1EA
	Final Drill D2.8-12mm for Ø3.3/4.1 Implant	3AA-190	1EA
	Final Drill D2.8-14mm for Ø3.3/4.1 Implant	3AA-191	1EA
	Final Drill D3.3-8mm for Ø4.8 Implant	3AA-192	1EA
	Final Drill D3.3-10mm for Ø4.8 Implant	3AA-193	1EA
8	Final Drill D3.3-12mm for Ø4.8 Implant	3AA-194	1EA
	Final Drill D3.3-14mm for Ø4.8 Implant	3AA-195	1EA
	Final Drill D4.0-8mm for Ø5.5 Implant	3AA-196	1EA
	Final Drill D4.0-10mm for Ø5.5 Implant	3AA-197	1EA
	Final Drill D4.0-12mm for Ø5.5 Implant	3AA-198	1EA
	Final Drill D4.0-14mm for Ø5.5 Implant	3AA-199	1EA
	Counter Sink D4.1	3AA-014	1EA
9	Counter Sink D4.8	3AA-015	1EA
	Counter Sink D5.5	3AA-016	1EA
	Profile Drill D3.5	3AA-065	1EA
10	Profile Drill D4.0	3AA-066	1EA
	Profile Drill D4.5	3AA-067	1EA
	Taps D3.3 for Ø3.3/4.1 Implant	3AA-017	1EA
	Taps D4.0 for Ø4.0/4.8 Implant	3AA-018	1EA
11	Taps D4.7 for Ø4.5/5.5 Implant	3AA-019	1EA
	Implant Driver Hex2.0 HP-L	3AA-056	1EA
	Implant Driver Hex2.0 RT-S	3AA-030	1EA
	Implant Driver Hex2.0 RT-L	3AA-039	1EA
12	Implant Driver Hex2.5 HP-L	3AA-057	1EA
	Implant Driver Hex2.5 RT-S	3AA-032	1EA
	Implant Driver Hex2.5 RT-L	3AA-040	1EA
	Screw Driver Hex1.25-HP-L	3AA-041	1EA
13	Screw Driver Hex1.25-RT-S	3AA-042	1EA
	Screw Driver Hex1.25-RT-L	3AA-043	1EA
14	Torque Ratchet 10-40Ncm	3AA-034	1EA


BIOMATE FULL KIT

- Color management:
- ◆ BIOMATE/ Blue(ø3.3) · Green(ø4.1) · Purple(ø4.8) · Orange(ø5.5)
 - ◆ BIOMATE PLUS/ Dark Blue(ø3.5) · Pink(ø4.0) · Brown(ø4.5)





Stopper for Drill D2.0/D2.8



Stopper for Drill D3.3/D4.0



Paralleling Pin




Torque Ratchet


	Description	Catalog No.	
1	BIOMATE FULL KIT	3AA-139	Full Instruments
2	Lance Drill D1.8	3AA-038	1EA
3	Initial Drill D2.0-L	3AA-007	1EA
4	Paralleling Pin D2.0- D2.8	3AA-052	3EA
5	Drill Extender	3AA-035	1EA
6	Final Drill D2.8-L for ϕ 3.3/4.1 Implant	3AA-009	1EA
	Final Drill D3.3-L for ϕ 4.8 Implant	3AA-011	1EA
	Final Drill D4.0-L for ϕ 5.5 Implant	3AA-013	1EA
9	Counter Sink D 4.1	3AA-014	1EA
	Counter Sink D 4.8	3AA-015	1EA
	Counter Sink D 5.5	3AA-016	1EA
10	Profile Drill D 3.5	3AA-065	1EA
	Profile Drill D 4.0	3AA-066	1EA
	Profile Drill D 4.5	3AA-067	1EA
11	Implant Driver Hex2.0 HP-S	3AA-029	1EA
	Implant Driver Hex2.0 HP-L	3AA-056	1EA
	Implant Driver Hex2.0 RT-S	3AA-030	1EA
	Implant Driver Hex2.0 RT-L	3AA-039	1EA
	Implant Driver Hex2.5 HP-S	3AA-031	1EA
	Implant Driver Hex2.5 HP-L	3AA-057	1EA
	Implant Driver Hex2.5 RT-S	3AA-032	1EA
12	Screw Driver Hex1.25-RT-S	3AA-042	1EA
	Screw Driver Hex1.25-RT-L	3AA-043	1EA
13	Stopper L8 for Drill D2.0/D2.8	3AA-021	1EA
	Stopper L10 for Drill D2.0/D2.8	3AA-022	1EA
	Stopper L12 for Drill D2.0/D2.8	3AA-023	1EA
	Stopper L14 for Drill D2.0/D2.8	3AA-024	1EA
	Stopper L8 for Drill D3.3/D4.0	3AA-025	1EA
	Stopper L10 for Drill D3.3/D4.0	3AA-026	1EA
	Stopper L12 for Drill D3.3/D4.0	3AA-027	1EA
14	Torque Ratchet 10-40Ncm	3AA-034	1EA

BIONavi MASTER KIT

Color management:

- ◆ BIOMATE/ Blue($\phi 3.3$) · Green($\phi 4.1$) · Purple($\phi 4.8$) · Orange($\phi 5.5$)
- ◆ BIOMATE PLUS/ Dark Blue($\phi 3.5$) · Pink($\phi 4.0$) · Brown($\phi 4.5$)





BIONAVI Implant Driver-HP



BIONAVI Implant Driver-RT



Screw Driver -HP

Screw Driver -RT



BIONAVI Taps D3.3 for Ø3.3/3.5/4.1 Implant

BIONAVI Taps D4.0 for Ø4.0/4.8 Implant


BIONAVI Taps D4.7 for Ø4.5/5.5 Implant



Torque Ratchet


	Description	Catalog No.
1	BIONAVI MASTER KIT	3AA-159 Full Instruments
2	BIONAVI Tissue Punch D3.0	3AA-N35 1EA
3	BIONAVI Bone Flattening	3AA-N03 1EA
4	BIONAVI Drill Tube D2.0/D5.3	3AA-N04 2EA
	BIONAVI Drill Tube D2.5/D5.3	3AA-N05 1EA
5	BIONAVI Lance Drill D2.0	3AA-N51 1EA
	BIONAVI Initial Drill D2.0-8mm	3AA-N24 1EA
6	BIONAVI Initial Drill D2.0-10mm	3AA-N25 1EA
	BIONAVI Initial Drill D2.0-12mm	3AA-N26 1EA
	BIONAVI Initial Drill D2.0-14mm	3AA-N27 1EA
	BIONAVI Final Drill D2.8-5mm for Ø3.3/3.5/4.1 Implant	3AA-N58 1EA
	BIONAVI Final Drill D2.8-8mm for Ø3.3/3.5/4.1 Implant	3AA-N06 1EA
	BIONAVI Final Drill D2.8-10mm for Ø3.3/3.5/4.1 Implant	3AA-N07 1EA
	BIONAVI Final Drill D2.8-12mm for Ø3.3/3.5/4.1 Implant	3AA-N08 1EA
	BIONAVI Final Drill D2.8-14mm for Ø3.3/3.5/4.1 Implant	3AA-N09 1EA
	BIONAVI Final Drill D3.3-5mm for Ø4.0/4.8 Implant	3AA-N59 1EA
	BIONAVI Final Drill D3.3-8mm for Ø4.0/4.8 Implant	3AA-N10 1EA
7	BIONAVI Final Drill D3.3-10mm for Ø4.0/4.8 Implant	3AA-N11 1EA
	BIONAVI Final Drill D3.3-12mm for Ø4.0/4.8 Implant	3AA-N12 1EA
	BIONAVI Final Drill D3.3-14mm for Ø4.0/4.8 Implant	3AA-N13 1EA
	BIONAVI Final Drill D4.0-5mm for Ø4.5/5.5 Implant	3AA-N60 1EA
	BIONAVI Final Drill D4.0-8mm for Ø4.5/5.5 Implant	3AA-N14 1EA
	BIONAVI Final Drill D4.0-10mm for Ø4.5/5.5 Implant	3AA-N15 1EA
	BIONAVI Final Drill D4.0-12mm for Ø4.5/5.5 Implant	3AA-N16 1EA
	BIONAVI Final Drill D4.0-14mm for Ø4.5/5.5 Implant	3AA-N17 1EA
	BIONAVI Counter Sink D4.1	3AA-N52 1EA
8	BIONAVI Counter Sink D4.8	3AA-N53 1EA
	Counter Sink D5.5	3AA-016 1EA
	BIONAVI Profile Drill D3.5	3AA-N30 1EA
9	BIONAVI Profile Drill D4.0	3AA-N31 1EA
	BIONAVI Profile Drill D4.5	3AA-N49 1EA
	BIONAVI Implant Driver-Stopper Hex 2.0-HP-S	3AA-N18 1EA
	BIONAVI Implant Driver-Stopper Hex 2.0-RT-S	3AA-N19 1EA
10	BIONAVI Implant Driver-Non Stopper Hex 2.0-RT	3AA-N22 1EA
	BIONAVI Implant Driver-Stopper Hex 2.5-HP-S	3AA-N20 1EA
	BIONAVI Implant Driver-Stopper Hex 2.5-RT-S	3AA-N21 1EA
	BIONAVI Implant Driver-Non Stopper Hex 2.5-RT	3AA-N23 1EA
11	BIONAVI Abutment Drill	3AA-N01 1EA
	Screw Driver Hex1.25-HP-L	3AA-041 1EA
12	Screw Driver Hex1.25-RT-S	3AA-042 1EA
	Screw Driver Hex1.25-RT-L	3AA-043 1EA
	BIONAVI Taps D3.3 for Ø3.3/4.1 Implant	3AA-N33 1EA
13	BIONAVI Taps D4.0 for Ø4.8 Implant	3AA-N34 1EA
	BIONAVI Taps D4.7 for Ø5.5 Implant	3AA-N50 1EA
14	Drill Extender-L	3AA-035 1EA
	BIONAVI Point Straight Drill-S	3AA-N28 1EA
15	BIONAVI Point Straight Drill-L	3AA-N29 1EA
16	Torque Ratchet 10-40Ncm	3AA-034 1EA

BIOMATE SL KIT



Color management:



- ◆ BIOMATE/ Blue(ø3.3) · Green(ø4.1) · Purple(ø4.8) · Orange(ø5.5)
- ◆ BIOMATE PLUS/ Dark Blue(ø3.5) · Pink(ø4.0) · Brown(ø4.5)



Description	Catalog No.	
1 BIOMATE SL KIT	3AA-063	Full Instruments
2 Lance Drill D1.8	3AA-038	1EA
3 Initial Drill D2.0-L	3AA-007	1EA
4 Paralleling Pin D2.0- D2.8	3AA-052	2EA
Final Drill D2.8-L for ø3.3/4.1 Implant	3AA-009	1EA
5 Final Drill D3.3-L for ø4.8 Implant	3AA-011	1EA
Final Drill D4.0-L for ø5.5 Implant	3AA-013	1EA
Counter Sink D4.1	3AA-014	1EA
6 Counter Sink D4.8	3AA-015	1EA
Counter Sink D5.5	3AA-016	1EA


Description	Catalog No.	
Profile Drill D3.5	3AA-065	1EA
7 Profile Drill D4.0	3AA-066	1EA
Profile Drill D4.5	3AA-067	1EA
Implant Driver Hex2.0 HP-L	3AA-056	1EA
8 Implant Driver Hex2.0 RT-L	3AA-039	1EA
Implant Driver Hex2.5 HP-L	3AA-057	1EA
Implant Driver Hex2.5 RT-L	3AA-040	1EA
9 Screw Driver Hex1.25-RT-L	3AA-043	1EA
10 Torque Ratchet 10~40Ncm	3AA-034	1EA

BIOMATE PROSTHETIC KIT




	Description	Catalog No.	
1	BIOMATE PROSTHETIC KIT	3AA-089	Full Instruments
2	Screw Driver Hex 1.25-HP-L	3AA-041	1EA
	Screw Driver Hex 1.25-RT-S	3AA-042	1EA
	Screw Driver Hex 1.25-RT-L	3AA-043	1EA
	Screw Driver Hex 1.25-RT-EL	3AA-148	1EA
3	Ball Abutment Driver Hex-HP	3AA-050	1EA
	Ball Abutment Driver Hex-RT-L	3AA-053	1EA
4	Torque Ratchet 10-40Ncm	3AA-034	1EA

POSITIONER / MULTI-UNIT KIT




Positioner Torque Driver



Multi-Unit Adapter



Screw Driver





Positioner Core Tool



Torque Ratchet (10-40Ncm)


	Description	Catalog No.	
1	Positioner Multi-Unit KIT	BSSIPM00000FA	Full Instruments
2	Positioner Torque Driver-S	3AA-085	1EA
	Positioner Torque Driver-L	3AA-086	1EA
3	Multi-Unit Abutment Adapter-HP	BSMUHP000000A	1EA
	Multi-Unit Abutment Adapter-RT	BSMURT000000A	1EA
4	Screw Driver Hex 1.25-RT-S	3AA-042	1EA
	Screw Driver Hex 1.25-RT-L	3AA-043	1EA
5	Positioner Core Tool	3AA-087	1EA
6	Torque Ratchet 10-40Ncm	3AA-034	1EA


BONE EXPANDER KIT



Description		Catalog No.	
1	BONE EXPANDER KIT	3AK-A00	Full Instruments
2	Crest Trimmer	3AK-A01	1EA
3	Saw	3AK-A02	1EA
4	Drill D1.8	3AK-A03	1EA
5	Expander D2.6	3AK-A04	1EA
	Expander D3.0	3AK-A05	1EA
	Expander D3.4	3AK-A06	1EA
	Expander D3.8	3AK-A07	1EA
	Expander D4.3	3AK-A08	1EA
6	Ratchet	3AK-A09	1EA
7	Ratchet Adapter	3AK-A10	1EA
8	Handpiece Adapter	3AK-A11	1EA
9	Hand / Straumann Adapter	3AK-A12	1EA

SINUS CRESTAL APPROACH KIT






Stopper 1mm-10mm



Bone Syringe

	Description	Catalog No.	
1	Sinus Crestal Approach Kit	3AK-B00	Full Instruments
2	Initial Drill D1.8	3AK-B13	1EA
	Initial Drill D2.3	3AK-B14	1EA
3	Sinus Reamer D2.8	3AK-B15	1EA
	Sinus Reamer D3.3	3AK-B16	1EA
	Sinus Reamer D3.8	3AK-B17	1EA
	Sinus Reamer D4.2	3AK-B18	1EA
4	Stopper 1mm	3AK-B01	1EA
	Stopper 2mm	3AK-B02	1EA
	Stopper 3mm	3AK-B03	1EA
	Stopper 4mm	3AK-B04	1EA
	Stopper 5mm	3AK-B05	1EA
	Stopper 6mm	3AK-B06	1EA
	Stopper 7mm	3AK-B07	1EA
	Stopper 8mm	3AK-B08	1EA
	Stopper 9mm	3AK-B09	1EA
	Stopper 10mm	3AK-B10	1EA
5	Aqua Lift Tap D2.8 (2.8~3.8)	3AK-B20	1EA
	Aqua Lift Tap D3.3 (3.3~4.2)	3AK-B21	1EA
6	Hand Adapter	3AK-B19	1EA
7	Diamond Bur D2.8	3AK-B11	1EA
8	Bone Syringe D3.5 / D4.0	3AK-B23	1EA
9	Depth Gauge	3AK-B22	1EA
10	Silicone Tube	3AK-B12	1EA

Biomate Plus SC Kit



Description	Dimension	Catalog No.
Biomate Plus SC Kit	Full Instruments 43PCS+1BOX	3AA-137

Biomate Master Kit



Description	Dimension	Catalog No.
Biomate Master Kit	Full Instruments 41PCS+1BOX	3AA-154

Biomate Full Kit



Description	Dimension	Catalog No.
Biomate Full Kit	Full Instruments 34PCS+1BOX	3AA-139

Biomate SL Kit









Description	Dimension	Catalog No.
Biomate SL Kit	Full Instruments 19PCS+1BOX	3AA-063

Biomate SS Kit




Description	Dimension	Catalog No.
Biomate SS Kit	Full Instruments 13PCS+1BOX	3AA-062

BIOnavi Master Kit	Description	Dimension	Catalog No.
	BIOnavi Master Kit	Full Instruments 51PCS+1BOX	3AA-159
BIOMATE Prosthetic Kit	Description	Dimension	Catalog No.
	Prosthetic Kit	Full Instruments 7PCS+1BOX	3AA-089
Positioner / Multi-Unit Kit	Description	Dimension	Catalog No.
	Positioner / Multi-Unit Kit	Full Instruments 8PCS+1BOX	BSSIPM00000FA
Sinus Crestal Approach Kit	Description	Dimension	Catalog No.
	Sinus Crestal Approach Kit	Full Instruments 23PCS+1BOX	3AK-B00
Bone Expander Kit	Description	Dimension	Catalog No.
	Bone Expander Kit	Full Instruments 12PCS+1BOX	3AK-A00
Implant & Screw Remover Kit (S)	Description	Dimension	Catalog No.
	Implant & Screw Remover Kit (s)	Full Instruments 24PCS+1BOX	3AK-D00


Surgical Instruments

Depth Marks on Biomate Instruments	74
Round Bur.....	75
Lance Drill.....	75
Lindemann Drill (Pilot Drill).....	75
Implant Driver	75
Initial Drill	76
Final Drill	76
Counter Sink	77
Profile Drill	77
Taps	77
Stopper For Drill	78
Drill Extender	78
Handpiece Adapter	78
Torque Ratchet.....	79
BioSmart Torque Ratchet	79
Paralleling Pin	79
Screw Driver	80
Depth Gauge	80
Ball Abutment Driver	80



Depth Marks on Biomate Instruments

Laser marks are made on the bladed end of the drills to indicate drilled depth for the practitioner




1. Initial Drill D2.0-L
2. Final Drill D2.0-L corresponds with Ø4.1mm/Ø3.5mm fixtures
3. Guide Pin D2.8mm
4. Taps D3.3mm corresponds with Ø 4.1mm/ Ø3.5mm fixtures
5. Ø 4.1mm fixtures with length of 8mm,10mm,12mm,14mm
6. Ø 3.3mm fixtures with length of 8mm,10mm,12mm,14mm
7. Profile Drill D3.5mm corresponds with Ø3.5mm fixtures

Surgical Instruments

Round Bur



Lance Drill




Name	TIP	Diameter(D)	Height(H)	Catalog No.
------	-----	-------------	-----------	-------------

Round Bur	—	D2.3	26	3AA-001
Lance Drill	0.5	D1.8	33	3AA-038

- Used in the initial stage of surgery to mark the position for implantation
- Used with rotation speed set at 1,200 rpm, torque 20 Ncm, feed water
- Round Bur can mark and smooth out a flat bone surface for drilling
- The pointed design on the Lance Drill provides a stable drilling into the cortical bone

Lindemann Drill (Pilot Drill)




Name	TIP	Diameter(D)	Height(H)	Catalog No.
------	-----	-------------	-----------	-------------

Lindemann Drill (Pilot Drill)	0.5	D2.0xL14	32	3AA-037
-------------------------------	-----	----------	----	---------


- Its side-cut design can correct the deviated position and angle of the initial drilling path
- Used with rotation speed set at 1,200 rpm, feed water
- Also called a Sidecut or a Lindemann

Implant Driver

D2.0



D2.5




Name	Diameter(D)	Height(H)	Catalog No.
------	-------------	-----------	-------------

Implant Driver	Hex2.0-HP-S	27	3AA-029
	Hex2.0-HP-L	32	3AA-056
	Hex2.0-HP-EL	37	3AA-124
	Hex2.0-RT-S	24	3AA-030
	Hex2.0-RT-L	26	3AA-039
	Hex2.0-RT-EL	29	3AA-126
	Hex2.5-HP-S	27	3AA-031
	Hex2.5-HP-L	32	3AA-057
	Hex2.5-HP-EL	37	3AA-125
	Hex2.5-RT-S	24	3AA-032
	Hex2.5-RT-L	26	3AA-040
	Hex2.5-RT-EL	29	3AA-127

- Hex D2.0 corresponds with D3.1/ 4.1 fixture to lock the fixture into the bone
- Hex D2.5 corresponds with D4.8/ 5.5 fixture to lock the fixture into the bone
- HP instruments are used with implant motor, RT instruments are used with the torque ratchet
- Line markings are 2mm apart to help determine the gingival height during minimal invasive surgery
- The dots indicate each flat surface of the internal hexagon of the fixture to help determine which direction the abutment is facing


Initial Drill



Name	TIP	Diameter(D)	Height(H)	Catalog No.
Initial Drill	0.5	D2.0-S	27	3AA-006
Initial Drill	0.5	D2.0-L	35	3AA-007

- The size of the initial drill is D2.0
- The measured length excludes the tip
- Used with rotation speed set at 1,200 rpm, torque 20 Ncm, feed water
- Used for drilling the hole to the required depth in the bone after locating point of implantation
- Initial Drill S is used when there is limited workspace in patient's mouth
- Initial Drill L can be used with a stopper when there are neighboring teeth


Final Drill







Name	TIP	Diameter(D)	Height(H)	Catalog No.	Fixture \varnothing
Final Drill	0.5	D2.8-S	27	3AA-008	$\varnothing 3.3$ $\varnothing 4.1$ $\varnothing 3.5$
Final Drill	0.5	D2.8-L	35	3AA-009	$\varnothing 3.3$ $\varnothing 4.1$ $\varnothing 3.5$
Final Drill	0.6	D3.3-S	27	3AA-010	$\varnothing 4.8$ $\varnothing 4.0$
Final Drill	0.6	D3.3-L	35	3AA-011	$\varnothing 4.8$ $\varnothing 4.0$
Final Drill	0.6	D4.0-S	27	3AA-012	$\varnothing 5.5$ $\varnothing 4.5$
Final Drill	0.6	D4.0-L	35	3AA-013	$\varnothing 5.5$ $\varnothing 4.5$


- Each drill corresponds to different diameters of fixtures
- The measured length excludes the tip
- Used with rotation speed set at 1,200 rpm, torque 20 Ncm, feed water
- Final Drill S is used when there is limited workspace in patient's mouth
- Final Drill L can be used with a stopper when there are neighboring teeth
- Final Drill D2.8 is used for enlarging the hole after using Initial Drill
- Final Drill D3.3 is used for enlarging the hole after using Final Drill D2.8
- Final Drill D4.0 is used for enlarging the hole after using Final Drill D3.3





Counter Sink




Name	Diameter(D)	Height(H)	Catalog No.	Fixture 
Counter Sink	D4.1	29	3AA-014	 Ø4.1
	D4.8	29	3AA-015	 Ø4.8
	D5.5	29	3AA-016	 Ø5.5
<ul style="list-style-type: none"> Used to trim the cortical bone with rotation speed 1,200rpm, torque 20Ncm, feed water; each drill corresponds to different diameters of fixtures Drill to the second laser mark for L10-L14 Implant or the first for L8mm Implant Counter Sink D4.1 is used after Final Drill D2.8 Counter Sink D4.8 is used after Final Drill D3.3 Counter Sink D5.5 is used after Final Drill D4.0 				









Profile Drill



Name	Diameter(D)	Height(H)	Catalog No.	Fixture 
Profile Drill	D3.5	28	3AA-065	 Ø3.5
	D4.0	28	3AA-066	 Ø4.0
	D4.5	28	3AA-067	 Ø4.5
<ul style="list-style-type: none"> Used to trim the cortical bone with rotation speed 1,200rpm, torque 20Ncm, feed water; each drill corresponds to different diameters of fixtures Profile Drill Profile Drill D3.5 is used after Final Drill D2.8 Profile Drill D4.0 is used after Final Drill D3.3 Profile Drill D4.5 is used after Final Drill D4.0 				



Taps




Name	Diameter(D)	Height(H)	Catalog No.	Fixture 
Taps	D3.3	29	3AA-017	 Ø3.3  Ø4.1  Ø3.5
	D4.0	29	3AA-018	 Ø4.8  Ø4.0
	D4.7	29	3AA-019	 Ø5.5  Ø4.5
<ul style="list-style-type: none"> Used on D1 bone to create threads inside the hole; each drill corresponds to different diameters of fixtures Used with rotation speed set at 20 rpm, torque 35 Ncm; set the implant motors to reversed rotation to withdraw the instrument after drilling Taps D3.3 is used after drilling with Final Drill D2.8 and Counter Sink D4.1 Taps D4.0 is used after drilling with Final Drill D3.3/ Profile Drill D4.0 Taps D4.7 is used after drilling with Final Drill D4.0/ Counter Sink D5.5 				

Stopper For Drill

D2.0 / D2.8




D3.3 / D4.0



Name	Diameter(D)	Height(H)	Catalog No.
Stopper For Drill L6	D2.0/D2.8	14	3AA-020
Stopper For Drill L8	D2.0/D2.8	12	3AA-021
Stopper For Drill L9	D2.0/D2.8	11	3AA-090
Stopper For Drill L10	D2.0/D2.8	10	3AA-022
Stopper For Drill L11	D2.0/D2.8	9	3AA-091
Stopper For Drill L12	D2.0/D2.8	8	3AA-023
Stopper For Drill L13	D2.0/D2.8	7	3AA-092
Stopper For Drill L14	D2.0/D2.8	6	3AA-024
Stopper For Drill L6	D3.3/D4.0	14	3AA-064
Stopper For Drill L8	D3.3/D4.0	12	3AA-025
Stopper For Drill L9	D3.3/D4.0	11	3AA-093
Stopper For Drill L10	D3.3/D4.0	10	3AA-026
Stopper For Drill L11	D3.3/D4.0	9	3AA-094
Stopper For Drill L12	D3.3/D4.0	8	3AA-027
Stopper For Drill L13	D3.3/D4.0	7	3AA-095
Stopper For Drill L14	D3.3/D4.0	6	3AA-028

- The Stopper is a safety sleeve that can be fit onto the Initial Drill-L or the Final Drill-L through the tip to prevent drilling too deep
- Stopper D2.0/2.8 is used with Initial Drill D2.0 and Final Drill D2.8
- Stopper D3.3/4.0 is used with Final Drill D3.3/4.0

Drill Extender



Name	Diameter(D)	Height(H)	Catalog No.
Drill Extender	L	27	3AA-035
	S	25	3AA-058

- Used for extending drills to avoid neighboring teeth

Handpiece Adapter




Name	Height(H)	Catalog No.
Handpiece Adapter	16	3AA-045

- Adapt instrument for use on handpiece to the torque ratchet
- Switching from motorized mode to manual mode
- How to use:




Torque Ratchet



Name	Diameter(D)	Height(H)	Catalog No.
Torque Ratchet	10-40Ncm	83	3AA-034

- Wrench to apply a constant torque (10/20/30Ncm) to screws and abutments
- When the set torque is applied, the neck of the Torque Wrench is bent for indication
- If a continuous force is applied while the neck is bent, excessive torque is applied, resulting in screw fracture
- Twist the adjustable end to set the required torque value; loosen fully to achieve infinite torque




BioSmart Torque Ratchet



Name	Diameter(D)	Catalog No.
BioSmart Torque Ratchet	80Ncm	BSSITR000000A

- A set of a two-way Torque Wrench and a Torque Connector
- Applying forward/reverse torque by rotating the Torque Wrench handle without removing the connector
- Applying torque according to the line marked with the torque value to be applied by pulling the bar
- Torque applied up to 80Ncm (15/30/45/80Ncm scale display)
- Washed and sterilized after use for storing


Paralleling Pin



Name	Diameter(D)	Height(H)	Catalog No.
Paralleling Pin	D2.0/2.8	18	3AA-052

- Placed in the hole to inspect the occlusion and the distance of the implant with neighboring teeth


Screw Driver



Name	Diameter(D)	Height(H)	Catalog No.
Screw Driver	Hex 1.25-HP-S	19	3AA-033
	Hex 1.25-HP-L	25	3AA-041
	Hex 1.25-HP-EL	31	3AA-145
	Hex 1.25-RT-SS	14.8	3AA-146
	Hex 1.25-RT-S	19.8	3AA-042
	Hex 1.25-RT-M	21.8	3AA-147
	Hex 1.25-RT-L	24.8	3AA-043
	Hex 1.25-RT-EL	31.8	3AA-148

- Used for connecting the healing components (Cover Screw/Healing Abutment) or the prosthetic components (Abutments) with the fixtures
- HP instruments are used with implant motor, RT instruments are used with the torque ratchet


Depth Gauge



Name		Height(H)	Catalog No.
Depth Gauge	—	139	3AA-044

- Used to inspect the height of the healing abutment and abutment after implant, and the depth of the hole

Ball Abutment Driver



Name	Diameter(D)	Height(H)	Catalog No.
Ball Abutment Driver	Hex-HP	18.5	3AA-050
	Hex-RT-S	20.0	3AA-051
	Hex-RT-L	26.0	3AA-053

- Used for Ball Abutment
- Hex-HP is used with handpiece at 30 Ncm
- Hex-HP is used with handpiece at 25 Ncm

Surgical Procedure

Preoperative Evaluation82

The Procedure.....83

Biomate Drilling Sequence of Instruments..... 89

Biomate Plus Drilling Sequence of Instruments91

Fixture Packaging and Label & Instrument Cleansing
and Maintenance93



*100% Insist
Biomate SWISS always do the best*

Preoperative Evaluation


Implant position of each case is different. Evaluation should be made in consideration of occlusion, neighboring teeth, thickness of the soft tissue, and the type of prosthetic components that will be used.

Distance between the bevel of the fixture and the neighboring teeth should be more than 1.5 mm

Distance between Biomate Implant and Neighboring Teeth				
Natural Teeth/Implant	3.3	4.1	4.8	5.5
Natural Teeth	3.2	3.6	3.9	4.3

Distance between Biomate Plus Implant and Neighboring Teeth			
Natural Teeth/Implant	3.5	4.0	4.5
Natural Teeth	3.3	3.5	3.8

Ø3.3 Prosthetic Platform




Osteotomy Center from Adjacent Tooth

Ø4.1 Prosthetic Platform




Osteotomy Center from Adjacent Tooth

Ø4.8 Prosthetic Platform




Osteotomy Center from Adjacent Tooth

Ø5.5 Prosthetic Platform




Osteotomy Center from Adjacent Tooth

If two fixtures are placed next to each other, the distance between the bevels of the two fixtures should be more than 3.0 mm



Measurement is Dependent on the Two Prosthetic Platform Diameters



Measurement is Dependent on the Two Prosthetic Platform Diameters

Distance between Biomate and Neighboring Teeth				
Fixture/ Fixture	3.3	4.1	4.8	5.5
3.3	6.3	6.7	7.1	7.4
4.1	6.7	7.1	7.5	7.8
4.8	7.1	7.5	7.8	8.2
5.5	7.4	7.8	8.2	8.5

Distance between Biomate-Plus Fixture and Neighboring Teeth			
Fixture/ Fixture	3.5	4.0	4.5
3.5	6.5	6.8	7
4.0	6.8	7	7.3
4.5	7	7.3	7.5

The Procedure

Operation Steps for Implanting Ø4.1 x L12 Fixture Operation Steps for Implanting Ø3.5 x L12 Fixture

STEP.1



Incision

Select a suitable scalpel to incise the gingiva and the periosteum at the desired implant site in order to expose the alveolar bone.

Note: The incision direction is subject to the patient's real bone condition, and the healing mode must also be considered (trans-gingival or submerged)

STEP.2



②

Lance Drill D1.8



Marking the Implant Position

After detaching mucosa, the Lance Drill is used to determine the implant site on the bone. The drilling depth is variable and is maximum the length of the implant.

Use the implant motor for the drilling procedure. Recommend speed max. 1,200 rpm (revolutions per minute) at 20Ncm; adjust appropriate water flow for cooling the drilling site.


Note: Surgical guide may be used to assist in marking the implant position.

STEP.3



⑥

Initial Drill D2.0(L)



Initial Drilling

Using the D2.8 mm Initial Drill, the implant length and axial alignment are then determined. The required drilling depth can be checked optically using the depth marking on the drills or using the optional depth pin. Recommended speed setting is 1,200 rpm at 20Ncm. The cavity is rinsed again with physiological saline solution.

Note: Stopper D2.0/2.8 may be used according to the needed depth.



STEP.4

STEP.5

STEP.6




⑦ ⑧ ⑨
Paralleling Pin
D2.0 / Ø8



⑪
Final Drill D2.8
TIP:0.5
8 10 12 14



⑬
Counter Sink D4.1
L10.L12.L14
L8



⑬
Profile Drill D3.5

Axial Alignment

Check the alignment with neighboring teeth and the occlusion by using paralleling pin. Alignment can be adjusted in subsequent steps.

Note: The incision direction is subject to the patient's real bone condition, and the healing mode must also be considered (trans-gingival or submerged)

Initial Extension Drilling

The initial drill hole is extended by using D2.8 mm final drill. Alignment can still be adjusted slightly at this stage. Recommended speed setting is 1,200 rpm at 20Ncm. Cooling is to use a chilled, sterile, physiological saline solution.

Note: Stopper D2.0/2.8 may be used according to the needed depth.

Expanding with Counter Sink (Only for BIOMATE)

Use Counter Sink D4.1 for trimming the cortical bone according to the patient's bone density. Enlarge the rim of the hole to correspond to the outer diameter of the fixture's platform.

Expanding with Profile Drill

Use Profile Drill D3.5 for trimming the cortical bone according to patient's bone density. Enlarge the rim of the hole to correspond to the outer diameter of the fixture's platform. Recommended speed setting is 1,200 rpm at 20 Ncm (feed water).

Note: Counter Sink is available for D1-D3 but not needed for D4 bone ; Profile Drill is available for D1-D2 bone , D3-D4 bone may be skipped with this step.

The Procedure

Operation Steps for Implanting Ø4.1 x L12 Fixture Operation Steps for Implanting Ø3.5 x L12 Fixture

STEP.7



STEP.8



Threading with Taps

BIOMATE fixture is applied with self-tapping design. However, for patients with high bone density (D1 bone), Taps D3.3 is Recommended. The reason of threading the drilled hole is to avoid excessive stress that might damage the bone and result a bone loss. Recommended speed is 20 rpm at 35 Ncm; set up a reversed rotation to withdraw tap after drilling.

Note: Taps is available for D1 bone, D2-D4 bone may be skipped with this step.

After drilling procedure, conduct implantation. Open the outer package to take out the sterile blister package. Peel the blister package open to obtain the fixture bottle and gently pull it open. Use Implant Driver D2.0-HP/RT to take the fixture out.



STEP:8



Note 1:

The hexagon and taper design of the fixture and the Implant Driver are made complimentary to each other. Gently press the Implant Driver to ensure it is firmly connected to the fixture before taking it out of the bottle.

Note 2:

Hold the Implant Driver with the fixture upside down to prevent the fixture from contacting other matter or dropping before placing it in the patient's mouth.

Note 3:

BIOMATE implant has a mount free design.



The Procedure

Operation Steps for Implanting Ø4.1 x L12 Fixture Operation Steps for Implanting Ø3.5 x L12 Fixture

STEP.9



STEP.10



Use the Implant Driver with implant motor or torque ratchet to screw the fixture into the bone with recommended torque 35Ncm. When the fixture cannot be fully screwed in, assess the necessity of unscrewing the fixture. Verify the diameter of the hole before a second approach.

Note: Using excessive force to screw in the fixture may damage the bone and cause bone infarction due to excessive stress. It is recommended to unscrew the fixture and re-drill the hole.




- Please turn Screw Driver D1.25-HP/RT(Tightened with 1.25 hex driver) counterclockwisely for 2 to 3 laps to take out the cover screw.

Cover Screw



STEP.11 Post-Implantation Procedure

Two Stage Surgery




- After implantation, use Screw Driver D1.25-HP/RT to take out the cover screw in the bottle cover. Ensure the Cover Screw is attached firmly with the Screw Driver to avoid the risk of dropping.
- Use the D1.25 screw driver to hand-tighten the Cover Screw into the fixture by torque ratchet(10Ncm).

Note 1: Do not use excessive force to prevent damaging the internal socket of the fixture.

Note 2: Healing abutment & cover screw, please keep the sterilization condition during surgery.

One Stage Surgery




- According to dentist's evaluation of patient's oral condition, a corresponding Healing Abutment can be placed right after implantation to omit incision a second time.
- Use the D1.25 screw driver to hand-tighten the Healing Abutment into the fixture by torque ratchet(10Ncm).


Note 1: Do not use excessive force to prevent damaging the internal socket of the fixture.

Note 2: Healing abutment & cover screw, please keep the sterilization condition during surgery.


Biomate Drilling Sequence of Instruments

D4	D3	D2	D1
A very thin layer of cortical bone surrounds a core of low-density trabecular bone Bone Quality	A thin layer of compact bone surrounds a core of dense trabecular bone of favorable strength	A thick layer of compact bone surrounds a core of dense trabecular bone	Almost the entire jaw is comprised of homogeneous compact bone
			
Soft Bone(D3~D4) Use a lower level Drill.	Medium Bone(D2~D3) *Profile Drill is recommended to remove cortical bone. *You can adjust the depth of Profile Drilling in accordance with bone density.	Hard Bone(D1~D2) *Tap Drill are recommended.	

Ø3.3 x L14 (Periodontal flap surgery)



Ø4.1 x L14 (Periodontal flap surgery)




• Recommended:


If you prefer to submerge the implant, we recommend Biomate Implant with 0.5-1.00 mm submerged.
If you follow up the drill mark, please exceed the mark on initial drill, final drill and counter sink in case of 0.5-1.00 mm submerged.

D4	D3	D2	D1
A very thin layer of cortical bone surrounds a core of low-density trabecular bone Bone Quality	A thin layer of compact bone surrounds a core of dense trabecular bone of favorable strength	A thick layer of compact bone surrounds a core of dense trabecular bone	Almost the entire jaw is comprised of homogeneous compact bone
Soft Bone(D3~D4) Use a lower level Drill.	Medium Bone(D2~D3) *Profile Drill is recommended to remove cortical bone. *You can adjust the depth of Profile Drilling in accordance with bone density.	Hard Bone(D1~D2) *Tap Drill are recommended.	

Ø4.8 x L14 (Periodontal flap surgery)




Ø5.5 x L14 (Periodontal flap surgery)




• Recommended:

If you prefer to submerge the implant, we recommend Biomate Implant with 0.5-1.00 mm submerged.
If you follow up the drill mark, please exceed the mark on initial drill, final drill and counter sink in case of 0.5-1.00 mm submerged.


Biomate Plus Drilling Sequence of Instruments




Ø3.5 x L10 (Periodontal flap surgery)



Ø4.0 x L10 (Periodontal flap surgery)




Ø4.5 x L10 (Periodontal flap surgery)




- **Recommended:** If you prefer to submerge the implant, we recommend Biomate Plus Implant with 0.5mm submerged. If you follow up the drill mark, please exceed the mark on initial drill, final drill and profile drill in case of 0.5mm submerged.

D4	D3	D2	D1
A very thin layer of cortical bone surrounds a core of low-density trabecular bone Bone Quality	A thin layer of compact bone surrounds a core of dense trabecular bone of favorable strength	A thick layer of compact bone surrounds a core of dense trabecular bone	Almost the entire jaw is comprised of homogeneous compact bone
Soft Bone(D3~D4) Use a lower level Drill.	Medium Bone(D2~D3) *Profile Drill is recommended to remove cortical bone. *You can adjust the depth of Profile Drilling in accordance with bone density.	Hard Bone(D1~D2) *Profile Drill and Tap Drill are recommended.	


※When the bone height is not enough to 8mm, it is recommended to use $\varnothing 4.5 \times L8$ implant, which makes the 1.0mm smooth surface appear above the bone plane. ($\varnothing 4.5 \times L7$)



$\varnothing 4.5 \times L8$ (Periodontal flap surgery)



$\varnothing 4.5 \times L8$ (Periodontal flap surgery)



- Recommended: If you prefer to submerge the implant, we recommend Biomate Plus Implant with 0.5mm submerged. If you follow up the drill mark, please exceed the mark on initial drill, final drill and profile drill in case of 0.5mm submerged.

Fixture Packaging and Label & Instrument Cleansing and Maintenance

Instrument Cleansing and Maintenance

Attention! Instruments are not sterilized when delivered, please autoclave prior to use.

- Please follow the instructions to clean and sterilize used surgical instruments.
01. During surgical process, soak used surgical instruments into saline solution.
- Cleaning** 02. After surgical process, use soft brush to clean remained blood stain and residues with clean water.
Soak surgical instruments into a container with quadruple enzyme cleaner that is covering all the instruments, and then put the container into a ultrasonic cleaning machine for 9 minutes. Finally, rinse the instruments with ultra-pure water several times and dry them.
-Do not use detergent containing aldehydes that can regulate and preserve protein.
-Please see manual of quadruple enzyme cleaner. (3M TM Rapid : Water, 1:100)
- Package** 03. Put the cleaned instruments back into surgical kit, and cover it with a surgical towel.
(Avoid collision among drills in order not to affect the cutting ability.)
- Sterilization** 04. Put the surgical kit with towel covered into sterilization pot (see manual of sterilization pot)
Recommended temperature: 132°C. Sterilization time: At least 4 minutes. Drying time: At least 30 minutes
- Storage** 05. After sterilization, keep the kit in a dust-proof and moisture-proof space. (Validity: Do not exceed 7 days)
- Notes:** Before using the surgical instrument, have to check. Discard the surgical instruments immediately if there is a defect as following:
-The blade becomes dull or damaged.
-Deformation (such as bending/twisting/folding)
-Surface corrosion

Recommendation for Use

To ensure quality use of instruments, instruments with cutting capability are recommended to be used less than 10 times. Please clean with sterile saline to prevent damage from excessive friction during surgery.

Sterilization Identification



Non Sterilized

Orange



Sterilized

Red

Fixture Package and Label


Fixture Product Packaging




Cover Screw



Implant internal label



Implant external label & sealing sticker



CATALOG

The Combination
of Technology
and Healthy Living



Biomate Swiss GmbH        **TFDA**

TEL: +886-7-951-1888 | 1F, NO.59, Luke 2nd Rd., Luzhu Dist., 821 Kaohsiung City, Taiwan
FAX: +886-7-727-8996 | www.biomateswiss.com info@biomateswiss.com

BI-IM-F025-A D21060101A